

Compact | Portable | Powerful

Copolymers have become increasingly commercially relevant due to the development of polymers with engineered properties. For example, acrylonitrile butadiene (ABS) for thermoplastics or styrene-butadiene (PS-PB) for car tires.

NMR Spectroscopy is a proven analytical technique for the characterization of fundamental polymer parameters because, in a single experiment you can:

- 1) Observe reaction completion & uniformity. Unreacted monomers are easily observable due to well-defined and resolved lines among the broad polymeric resonances.
- 2) Quantify the relative composition through pre-defined integral regions.
- 3) Quantify the relative percentages of structural isomers (e.g., branched vs. linear).
- 4) Observe tacticity and stereoisomers
- 5) Determine molecular weight, molecular number and polydispersity index (PDI)

Polymer characterization has been performed on traditional NMR systems for years, but now this technique is also available directly the benchtop. The affordable, lightweight and intuitive NMReady-60 is perfectly suited to improve analytical workflow. This high-resolution 60 MHz spectrometer generates results that are comparable and statistically equivalent to a high field super-conducting spectrometer. Unlike traditional spectrometers, however, the NMReady-60 is low maintenance, robust and easy-to-use.

## NMReady-60 Benchtop NMR Spectrometer for Block Copolymer Characterization



### SPECIFICATIONS:

Frequency:  $^1\text{H}$  60 MHz,  $^{31}\text{P}$  24 MHz  
Magnet: permanent, cryogen-free  
Stray field: confined within enclosure  
Uses Standard 5 mm NMR Tubes  
Spectral resolution:  
LW at 50%: 1.5 Hz  
LW at 1%: <25 Hz  
Power: 100-240 VAC, 50-60 Hz  
Weight: 50lbs  
Dimensions: 9.5 x 11 x 17"  
User Interface: Built-in Touchscreen  
Connectivity: USB, Ethernet

[www.nanalysis.com](http://www.nanalysis.com)



