


Determining Isotopic Ratio by Benchtop NMR: ^{10}B vs ^{11}B



In NMR, some nuclei have more than one NMR active isotope. ^{11}B (boron-11) and ^{10}B (boron-10) are the two stable isotopes of boron which are both NMR active. The ^1H NMR spectrum of sodium borohydride can be used indirectly to observe both nuclei in a single spectrum through their coupling to proton. The relative percentages of each isotope can be calculated from integrations. The spectrum shows a quartet (for $\text{Na}^{11}\text{BH}_4$) and septet (for $\text{Na}^{10}\text{BH}_4$). In the presence of a boron atom, more electropositive compared to a typical methyl carbon, the protons are shielded and appear at a lower frequency at about $\delta = 0.03$ ppm in D_2O . The spectrum shows that the experimentally determined ratio of 81.2 : 18.8 corresponds quite well to the known ratio of 80.1 : 19.9. Determining isotopic ratios of different nuclei demonstrates one of the uses of NMR beyond structural characterization.

$$2nI + 1$$

Splitting and multiplicity rule
where n is the number of coupled nuclei, and I is their nuclear spin




Sodium borohydride

^1H NMR
(60 MHz, D_2O)


^{10}B
septet
1:1:1:1:1:1:1
($^1J_{\text{BH}} = 27.0$ Hz)

^{11}B
quartet
1:1:1:1
($^1J_{\text{BH}} = 80.6$ Hz)

Nuclei	Nuclear Spin, I	Natural Abundance	NMR Active
^{10}B	3	19.9%	✓
^{11}B	$3/2$	80.1%	✓



- Protons
- Neutrons
- Electrons



Quadrupolar nuclei
have a nuclear spin greater than $1/2$

Read more about the isotopic ratio of boron through benchtop NMR here:

<https://www.nanalysis.com/nmready-blog/2021/11/2/eat-your-heart-out-mass-spec-measuring-10b11b-isotopic-ratio-by-nmr-spectroscopy>

https://www.quantanalitica.com/wp-content/uploads/2021/06/NMR_B-isotopes_ENG.pdf