



## BIOCHEMISTRY

### UNDERGRADUATE EXPERIMENT

# Determination of Reaction Kinetics for Hydrolysis of *N*-acetyl-DL-methionine

[www.nanalysis.com](http://www.nanalysis.com)

## INTRODUCTION

Enzymes are proteins that catalyze the organic reactions vital to sustain living organisms. The enzymatic reaction begins when the substrate (**S**) reversibly binds to the active site of the enzyme (**E**) to form an enzyme-substrate complex (**E-S**) with rate constants of  $k_1$  and  $k_{-1}$ . This is followed by the second step where the enzyme releases the product (**P**) with a rate constant of  $k_2$ .<sup>[1]</sup> The general reaction scheme of an enzyme catalyzed reaction is shown below.



To further understand the behavior of enzymes, a kinetic description of their activity is essential. One of the best-known models of enzyme kinetics is the Michaelis-Menten model.<sup>[2]</sup> The model is defined by an equation that relates the reaction rate,  $v$  (i.e. the rate of the formation of [**P**]), to the concentration of the substrate, [**S**]. The Michaelis-Menten equation is given below:

$$v = \frac{d[P]}{dt} = \frac{V_{max}[S]}{K_M + [S]}$$

From the Michaelis-Menten model, two important parameters can be determined,  $V_{max}$  and  $K_M$ .  $V_{max}$  represents the maximum rate of product formation at a saturating substrate concentration and is a measure of the efficiency of the enzyme as a catalyst. The Michaelis constant,  $K_M$ , represents the concentration of substrate at which the reaction rate is half of  $V_{max}$  and is often used to quantify the affinity of the active site for the substrate (the smaller the  $K_M$  value the higher the affinity). Typically,  $V_{max}$  and  $K_M$  are obtained by determining the initial reaction rate of an enzyme at varying substrate concentrations.<sup>[3]</sup> The reaction rate is then plotted against concentration to generate a Michaelis-Menten plot. By reciprocating both axes on the Michaelis-Menten plot, the Lineweaver-Burk plot can be obtained from which the  $V_{max}$  and  $K_M$  can be extracted from the line of best fit.

In this experiment, adapted from a J. Chem. Ed. article published by Olsen and Giles,<sup>[4]</sup> the enzymatic hydrolysis of *N*-acetyl-L-methionine by porcine acylase (*N*-acyl-L-aminoacid amidohydrolase) is studied. This reaction can be readily monitored via <sup>1</sup>H NMR spectroscopy with the NMReady-60. The data obtained from a single reaction can then be used to construct both a Michaelis-Menten and Lineweaver-Burk plot for a fast and semi-quantitative enzyme kinetics analysis.

## PROCEDURE


### Preparing Stock Solutions

*N*-acetyl-DL-methionine (0.382 g) was suspended in 2 mL of D<sub>2</sub>O along with 0.112 g of KH<sub>2</sub>PO<sub>4</sub>. Sodium hydroxide (2 M solution in D<sub>2</sub>O) was added carefully to bring the pH to 7 using pH paper. The resulting solution is then diluted to 5 mL in a volumetric flask using D<sub>2</sub>O. The final solution contained 400 mM of *N*-acetyl-DL-methionine. A stock solution of enzyme is prepared by dissolving 10 mg of porcine acylase and 1.5 mg of CoCl<sub>2</sub>·6H<sub>2</sub>O in 10 mL of D<sub>2</sub>O.


### Monitoring the Reaction with the NMReady-60

The solution of *N*-acetyl-DL-methionine (500 μL) is transferred to an NMR tube and a <sup>1</sup>H NMR spectrum was obtained (spectral width = 20 ppm, spectral centre = 5 ppm, number of scans = 16, delay = 0.5 sec, number of points = 4096). The reaction is initiated by adding 100 μL of the enzyme solution to the NMR tube followed by vigorous mixing. A <sup>1</sup>H NMR spectrum is recorded every 4 minutes for 2 hours using the kinetics module on the NMReady-60 (wait type = linear, number of clusters = 40, wait units = seconds, wait time (tau) = 160). To monitor the progress of the reaction, the integrals of the α-methine protons were measured for the reactant (*N*-acetyl-DL-methionine, 4.25 ppm) and product (L-methionine, 3.85 ppm).

## RESULTS




**Figure 1.** Stacked plot of  $^1\text{H}$  NMR spectra of the hydrolysis of *N*-acetyl-DL-methionine by porcine acylase to produce L-methionine.



**Figure 2.** Plot of substrate concentration over time of the reaction.

**Figure 3.** Michaelis-Menten plot of the reaction. The data was fitted to:  
$$V = (V_{\max} [S]) / (K_M + [S])$$



**Figure 4.** Lineweaver-Burk plot of the reaction. The data was fitted to the equation  $1/V = (K_M/V_{\max} [S]) + 1/V_{\max}$  from which the values of  $K_M$  (0.24 mol L $^{-1}$ ) and  $V_{\max}$  (0.3152 mmol L $^{-1}$  s $^{-1}$ ) were extracted.

## DISCUSSION

As seen in Figure 1, the  $^1\text{H}$  NMR spectrum of the hydrolysis reaction shows the depletion of the substrate, *N*-acetyl-DL-methionine (4.25 ppm), and the simultaneous appearance of the product, L-methionine (3.85 ppm). It is seen that the signal at 4.25 ppm never completely disappears because the D-enantiomer of the racemic mixture remains in the solution and does not get hydrolyzed by the porcine acylase. Figure 2 displays the plot of substrate concentration over time. The reaction is complete within an hour as the substrate concentration reaches a plateau. In Figure 3, the Michaelis-Menten plot illustrates the change of reaction rate as a function of substrate concentration. While the Michaelis-Menten experiment is typically carried out by measuring the reaction rate at several initial substrate concentrations, the experiment is condensed into one reaction in this case. By acquiring multiple  $^1\text{H}$  NMR spectra as the reaction proceeds, the substrate concentration can be determined from each spectrum and the reaction rate can be approximated by calculating the change in substrate concentration over a known time interval. Therefore, at higher substrate concentration it is seen that the reaction rate begins to reach a plateau which represents the  $V_{\max}$  at this substrate concentration. From the Michaelis-Menten plot, the Lineweaver-Burk plot (Figure 4) is constructed by reciprocating both axes. Subsequently, it was found that the  $K_M = 0.24$  mol L $^{-1}$  and  $V_{\max} = 0.3152$  mmol L $^{-1}$  s $^{-1}$ .

## CONCLUSIONS

In this experiment the enzymatic hydrolysis of *N*-acetyl-L-methionine was studied. Due to the difference in chemical shifts of the  $\alpha$ -methine protons in the substrate and product,  $^1\text{H}$  NMR spectroscopy could be used to monitor the progress of the reaction using the NMReady-60 instrument. Furthermore, quantitative data was obtained from the spectra that was used to construct a Michaelis-Menten and Lineweaver-Burk plot which were then used to determine the  $V_{\max}$  and  $K_M$  values of the enzymatic reaction.

## REFERENCES

- [1] Le, H.; Algaze, S.; Tan, E. Michaelis-Menten Kinetics [https://chem.libretexts.org/Textbook\\_Maps/Biological\\_Chemistry/Catalysts/Enzymatic\\_Kinetics/Michaelis-Menten\\_Kinetics](https://chem.libretexts.org/Textbook_Maps/Biological_Chemistry/Catalysts/Enzymatic_Kinetics/Michaelis-Menten_Kinetics) (accessed Dec 4, 2018).
- [2] Blanco, A.; Blanco, G. Medical biochemistry; Academic Press: London, United Kingdom, **2017**; pp. 153-175.
- [3] Berg, J.; Tymoczko, J.; Stryer, L. Biochemistry; 5th ed.; W.H. Freeman and Co.: New York, **2002**.
- [4] Olsen, R., Olsen, J. and Giles, G. "An Enzyme Kinetics Experiment for the Undergraduate Organic Chemistry Laboratory." *J. Chem, Educ.*, **2010**, 87(9), pp.956-957.

## DATA ACCESSIBILITY

The data can be processed directly on the NMReady-60 and printed and/or exported directly to a USB or networked file where it can be worked up using third party NMR processing software.

For additional ideas of how to incorporate the NMReady-60<sup>TM</sup> benchtop NMR spectrometer into undergraduate laboratories please see:

- 1) pH, pK<sub>a</sub> and Chemical Shift
- 2) Isomerization of Mo complexes via  $^{31}\text{P}$  NMR Spectroscopy
- 3) Aldol Condensation

or visit

[www.nanalysis.com/sample-experiments](http://www.nanalysis.com/sample-experiments)

**CONNECT WITH US**

 [www.nanalysis.com](http://www.nanalysis.com)

 1.855.NMReady

 @nanalysis

